


APP-CHEM NEWSLETTER



App-Chem will be showing at the **AquaTech Amsterdam** Exhibition; our stand is situated in The Analytical Hall. The show will see us re-launching a number of old favourites and a range of new and innovative products. We will also be looking to expand

our distribution network throughout the world enabling us to supply on a global basis.

Please Note:

As part of App-Chem's commitment to continued investment the website is going through a transition period which will be completed by the New Year.



Burette Kit

The only kit on the market to take 5 Burettes

This kit offers the user the ultimate testing facility in the field.

The use of the burettes allows for:

- Better accuracy
- Simple test identification
- Speedy set up
- Easy cleaning



The case you want is back!

App-Chem re-launch their Town Test Kit with **drop-down work surface** as standard. This rugged case will be 'book opening' (see photo) with lockable hinges. The case can be customised to your requirements as per all App-Chem Kits.



Titration



- New electronic rechargeable titrator
- Multitude of functions
- Positive displacement
- Use with standard cartridge
- Accurate multiple sample dispensing
- Can take cartridges from 0.5-50ml
- Large backlit LCD

ALTERNATIVELY APP-CHEM'S **MANUAL TITRATION SYSTEM** IS AVAILABLE. USABLE OVER THE ENTIRE RANGE BUT WITHOUT THE ADDITIONAL COST OF A TITRATOR. PLEASE ASK FOR MORE DETAILS.



HQ40d-input multi-parameter

Plug-and-play with any two IntelliCAL probes to measure pH, conductivity or dissolved oxygen using LDO technology in any combination. The HQ40d meter has three connection ports - two ports for simultaneous measurement using two IntelliCAL probes and an additional port for USB and external power.



Measurements can be made and viewed simultaneously. Transfer the data log or individual readings (GLP/ISO) to PC or USB memory stick. The meter is password-protected for accurate, tamper-proof reporting. The HQ40d meter comes with a field kit including carry case, protective gloves and sample cups.

- **Ready to use out of the box**
- **Calibration History stored in the probe itself**
 - Move probes between meters without the need to re-calibrate.
- **Versatile and durable**
 - All IntelliCAL probes are fully submersible for up to 24 hours.
- **Two probe styles to choose from**
 1. IntelliCAL Standard Probes
 2. IntelliCAL Rugged Probes

Getting in touch

Whether you want to place an order, request a catalogue, have a query, or would like to arrange a visit please get in touch with us in one of the following ways.

graham.allen@app-chem.com
catherine.bennion@app-chem.com
aimee.williams@app-chem.com
wendy.evans@app-chem.com
susan.bentley@app-chem.com

UK OFFICE

App-Chem House
Hugmore Lane
Llan y Pwll
Wrexham
LL13 9YE



TELEPHONE/ FAX

TELEPHONE: +44 (0) 1978 661551
FAX: +44 (0) 1978 661611

Managing Director
Technical Sales Manager
Sales Administrator
Sales Administrator
Accounts Manager

www.app-chem.com